




Data Sheet

Anti-Human-HLA-B,C

Catalog Number : MA010 Purified-anti- human HLA-**B,C**
MA010.1 Biotin- anti- human HLA-**B,C**
Specificity: Human beta-2-microglobulin
Buffer Solution: PBS pH 7.2, 1 mg/ml BSA, 0.09% NaN₃.
Storage Conditions: Store at 4°C. Avoid repeated freeze/thaw cycles.
Applications Tested



45kD

Balb/C-A2

Balb/C-B27

HLA- A2 transgene mouse spleen cell (left) and HLA-B27 transgene mouse spleen cell (right) were lysed and separated on SDS-PAGE gel. Anti-HLA-B,C was used to detect the existence of HLA B27 heavy chain. This antibody do not react with HLA-A, it react with HLA B and weak react with HLA-C.

Related Products

Cat	Description	Label
SA001	Anti-Mouse Ig	HRP
SA002	Anti-Rabbit Ig	HRP
SA003	Anti-Goat Ig	HRP
PD011	Streptavidin	HRP

For research purposes only. Not for therapeutic or diagnostic use.

www.epigenbiotec.com

Tel:010-64807408, 013910899152

epigen.info@gmail.com